

# STIMULATING NEWS

## AMERICAN IMEX

OFFERING THE FINEST PRODUCTS AT COMPETITIVE PRICES FOR OVER 25 YEARS  
QUARTERLY NEWSLETTER; FALL 2007

### YOUR INFORMATION PORTAL FOR ELECTROTHERAPY



[WWW.AMERICANIMEX.COM](http://WWW.AMERICANIMEX.COM)

- Full Product Descriptions
- Product Specification Sheet downloads
- Electrodes & Conductive Garments
- **Placement Guides**
- **Reimbursement Codes & Information**
- Product Specials and Promotions

### Sterile Electrodes

Post-operative use of TENS and Interferential is on the rise. Sterile Surgical Electrodes are a must for placement under a surgical dressing. We now manufacture the #40010-535 Sterile 2" round electrode. The 4 pack of electrodes with 9" color-coded (red, black) lead wires for easier access after application of the dressing.

Dealer Price \$153.00\* per case of 10 packs

\*1-7 case price

American Imex • 16520 Aston Street • Irvine CA 92606 • Tel: (800) 521-8286 • Fax: (949) 852-1245

### NEW A-TENS+

#### + MEANS MORE FEATURES

The New & Improved A-TENS+ is now available.



As buyers of the original A-TENS have come to expect; the A-TENS+ continues to offer the features that made it great; with improved circuitry, improved battery life and a fresh new look. Full Function

TENS- Conventional, Burst, Modulation Width and Modulation Rate PLUS 2 new Dense Disperse Modes for maximum effectiveness.

### INTERFERENTIAL CURRENT

The Interspec IF II combines dual channel Interferential with an EMS mode for Post Operative pain and spasm control. Six pre-set programs and one fully customizable program allow quick selection of the protocol that best



fit the patient's needs. Features include sequential IF/EMS stimulation and Vectoring IF. Stimulation time can be set from 1-30 minutes or continuously for severe post-surgical joint and deep tissue pain.

Full digital interface allows return to all settings used in previous sessions. AC adaptor included.

**CALL NOW!**

**(800) 521-8286**