

# MicroCare II

MicroCurrent Stimulator

Switchable from subsensory microcurrent stimulation to traditional milliampere conventional TENS stimulation. Choose from two waveforms and four prest modulation frequencies for maximum results.

Channel:	Dual, isolated between channels
Carrier Frequency**	4,000 Hz
Output:	High: 0-9 Milliamperes Low: 0-900 Microamperes
Modulated Frequency:	High: 292 Hz Moderate: 30 Hz Low: 0.3 Hz
Power Source:	9-Volt Battery.
Dimensions:	3.13" (93 mm) x 2.5" (64 mm) x 1"
Weight:	4.25 oz (119 grams (battery included)).



## **AMERICAN IMEX**

16520 Aston Street, Irvine CA 92606 USA  
949-553-8885 800-521-8286  
fax: 949-852-1245  
www.americanimex.com info@americanimex.com