


NuTrac™ Pelvator

The NuTrac Pelvator is intended to provide electrical stimulation and neuromuscular re-education for the purpose of rehabilitation of weak pelvic floor muscles for the treatment of stress, urge and mixed urinary incontinence in women

The NuTrac™ Pelvator is an advanced and modern digital dual channel STIM unit featuring 11 built-in continence programs and 3 custom programs. The NuTrac™ Pelvator unit will help with a wide variety of continence problems and is widely used by therapists around the world.


PELVATOR

- 11 Built-in programs;
- 3 Customizable, 5 segment programs
- LCD displays Intensity, Rate, Pulse Width, Time and program Mode
- Time setting up to 90 minutes
- Constant Current (extends battery life)
- Work / Rest periods from 2 - 99 seconds
- Range: Rate: 2 Hz - 100 Hz, Pulse Width: 50 μ S to 450 μ S
- Lock mode function for recording time in use and locking the units parameters
- Automatic switch off after 4 minutes if unit is not in use
- Suited for the experienced therapist while still maintaining its ease of use



VAGINAL PROBE

- Weight:42 grams
- Dimensions:length 87 mm,diameter 26 mm
- Body of probe:biocompatible ABS plastic



AMERICAN IMEX